

## REHABILITATION FOOTWEAR

### MEDICAL-SURGICAL SHOE

- Patented metatarsal shank provides control of forefoot motion to protect surgery to the metatarsal head area of the forefoot.
- Reinforced padded heel grips the foot reducing slippage



	Women		Men		Pediatric	
Size	Reference	Size	Reference	Size	Reference	Size
Small	DAR-MSW1N	4 - 6	DAR-MSM1N	6 - 8	DAR-MSPN	12 - 1
Medium	DAR-MSW2N	6 1/2 - 8	DAR-MSM2N	8 1/2 - 10		
Large	DAR-MSW3N	8 1/2 - 10	DAR-MSM3N	10 1/2 - 12		
X-Large			DAR-MSM4N	12 1/2 - 14		

### MEDICAL-SURGICAL SHOE WITH SQUARE TOE

- Rocker sole reduces plantar pressure
- Zoned outsole provides traction under the mid-foot
- MetaShank™ protection provides more rigid control under the metatarsal heads.
- Unique forefoot closure provides uniform security without the pressure points of buckles and straps
- The ankle strap seats the foot firmly in the sole and eliminates heel friction
- Removable insole fully customizable with over twice the padding
- The square toe design provides additional protection and better left/right fit



	Women		Men		Pediatric	
Size	Reference	Size	Reference	Size	Reference	Size
Small	DAR-MQW1B	4 - 6	DAR-MQM1B	6 - 8	DAR- MQPB	12 - 1
Medium	DAR-MQW2B	6 1/2 - 8	DAR-MQM2B	8 1/2 - 10		
Large	DAR-MQW3B	8 1/2 - 10	DAR-MQM3B	10 1/2 - 12		
X-Large			DAR-MQM4B	12 1/2 - 14		

### PEG-ASSIST™ INSOLE

- Inexpensive, time saving product to effectively off-load plantar pressures in the medical-surgical shoe with square toe.
- Made of a combination of Plastazote®, Poron® and EVA.
- Available in 7 sizes to fit the medical-surgical shoe with square toe above (DAR-MQ Series).

Size	Women	Men
Small	DAR-PTQW1	DAR-PTQM1
Medium	DAR-PTQW2	DAR-PTQM2
Large	DAR-PTQW3	DAR-PTQM3
X-Large		DAR-PTQM4




## REHABILITATION FOOTWEAR

### ORTHOWEDGE™ SHOE

Clinically proven to effectively reduce weight bearing pressure from the forefoot to promote proper healing.

	Shoe Size		
Size	Reference	Women	Men
X-Small	DAR-OW0B	4 - 7	
Small	DAR-OW1B	7 1/2 - 10	6 - 8
Medium	DAR-OW2B	10 1/2 - 13	8 1/2 - 10
Large	DAR-OW3B		10 1/2 - 12
X-Large	DAR-OW4B		12 1/2 - 14



### HEELWEDGE™ SHOE

- Off-loads the heel by shifting weight to the mid-foot to promote faster healing after surgery, trauma or when wounds are present on the heel.
- The ankle strap seats the foot firmly in the sole and eliminates heel friction.
- The square toe design provides additional protection and better left/right fit.

	Shoe Size		
Size	Reference	Women	Men
X-Small	DAR-HW0B	4 - 7	
Small	DAR-HW1B	7 1/2 - 10	6 - 8
Medium	DAR-HW2B	10 1/2 - 13	8 1/2 - 10
Large	DAR-HW3B		10 1/2 - 12
X-Large	DAR-HW4B		12 1/2 - 14

### WOUND CARE SHOE SYSTEM

The ultimate for off-loading lesions of the foot to permit healing of diabetic foot ulcers and infections.

Rocker sole to provide easy ambulation

- Easy hook and loop closure
- Sold in pairs to avoid sole height problems with the unaffected foot
- Allow health provider to off-load areas of the foot while the patient remains ambulatory
- Permits offloading ulcers, infected areas, bony prominences on both the plantar surface of the foot and protection of other areas of the foot



	Shoe Size		
Size	Reference	Women	Men
X-Small	DAR-WCSS0B	4 - 5 1/2	
Small	DAR-WCSS1B	6 - 8 1/2	6 - 7 1/2
Small/Med	DAR-WCSS2B	9 - 9 1/2	8 - 8 1/2
Medium	DAR-WCSS3B	10 - 10 1/2	9 - 9 1/2
Large	DAR-WCSS4B	11 - 12 1/2	10 - 11 1/2
X-Large	DAR-WCSS5B	13 - 14 1/2	12 - 13 1/2
XX-Large	DAR-WCSS7B	15 - 17	14 - 16

## REHABILITATION FOOTWEAR

### ALL PURPOSE BOOT

- Can be used as a post operative shoe and a cast boot
- The removable EVA insole is fully customizable to offload pressure points if needed.
- A wider sole and deep opening gives more room to accommodate large dressings, bandages and casts
- A wide ankle strap eliminates heel slippage as well as seats and secures the foot or cast in the shoe
- The rocker bottom encourages smooth roll through gait pattern and reduces plantar pressure .
- MetaShank™ construction provides for a reinforced sole under the metatarsal heads to reduce motion and the strapless forefoot closure provides uniform security without the pressure of buckles or straps



Reference	Size	Shoe Size	
		Women	Men
DAR-APB-XS	X-Small	4 - 7	-
DAR-APB-SM	Small	7 1/2 - 10	6 - 8
DAR-APB-MD	Medium	10 1/2 - 13	8 1/2 - 10
DAR-APB-LG	Large	-	10 1/2 - 12
DAR-APB-XL	X-Large	-	12 1/2 - 14

### STANDARD OPEN TOE CAST BOOT

Economical canvas cast boot with 3 hook and loop closures. Rocker bottom of 3 layers of EVA foam provides shock absorption. Canvas upper.

Reference	Size	Shoe Size	
		Women	Men
DAR-CBL00	X-Small	0 - 5	-
DAR-CBL01	Small	5 1/2 - 7 1/2	0 - 6
DAR-CBL02	Medium	8 - 10	6 1/2 - 8 1/2
DAR-CBL03	Large	-	9 - 11
DAR-CBL04	X-Large	-	12 - 14





## REHABILITATION FOOTWEAR

### BODY ARMOR™ CAST SHOE

- Soft but durable EVA provides shock absorption while protecting the cast from the elements.
- Unique bungee closure system is easy to use, easy to adjust and prevents the heel from slipping out of the cast shoe.
- The low profile, modern design will appeal to most patients.



Reference	Size	Shoe Size	
		Women	Men
DAR-BACS1	X-Small	3 - 5	-
DAR-BACS2	Small	5 1/2 - 8	-
DAR-BACS3	Medium	8 1/2 - 10	6 - 8
DAR-BACS4	Large	10 1/2 - 13	8 1/2 - 11
DAR-BACS5	X-Large	13 +	11.5 - 14

### BODY ARMOR™ WALKER

- Cast replacement walker provides rigid contoured upper with removable vented liner to provide protection, immobilization, compression for foot and ankle conditions where ambulation is permitted.
- Adjustable "ski boot" closure.

Reference		Size	Shoe Size	
Short	Standard		Women	Men
DAR-BAWS-SM	DAR-BAWH-SM	Small	5 1/2 - 8	-
DAR-BAWS-MD	DAR-BAWH-MD	Medium	8 1/2 - 10	6 - 8
DAR-BAWS-LG	DAR-BAWH-LG	Large	10 1/2 - 13	8 1/2 - 11



### BODY ARMOR™ WALKER LINER

Reference	Size
DAR-BAWLIN-SM	Small
DAR-BAWLIN-MD	Medium
DAR-BAWLIN-LG	Large

## REHABILITATION FOOTWEAR

### BODY ARMOR™ WALKER II

- Second generation walker with built in vented liner and raised cooling vents for greater comfort.
- Removable multi-density insoles that allow offloading with easy on and off for cleansing, wound care and dressing change.
- Flexible, upper conforms to leg to compress the foot and ankle.
- Locking ratcheted closures and raised air vents allow the escape of body heat.
- Adjustable circumferential compression reduces edema and promotes healing.
- Available in standard and short lengths.



Reference		Shoe Size		
Standard Length	Short Length	Size	Women	Men
DAR-BAT0	DAR-BATS0	X-Small	0 - 4	0 - 3
DAR-BAT1	DAR-BATS1	Small	4 1/2 - 6 1/2	3 1/2 - 6
DAR-BAT2	DAR-BATS2	Medium	7 - 10	6 1/2 - 9
DAR-BAT3	DAR-BATS3	Large	10 1/2 - 13	9 1/2 - 13
DAR-BAT4	DAR-BATS4	X-Large	-	13 1/2 - 17