

## MobiDig® Flexion Finger Splint

Finger splint for flexion of DIP or PIP joints.  
By careful bending of spring, flexion can be increased or decreased.  
High wearing comfort by elaborated, wear-resistant, foam padding.  
Linear spring force.


### Indications

Flexion of PIP (L1) or DIP (L2) joints

### Sizes

Widths as shown ( 1cm = 0.3937 in).

Reference	Size	L1 or L2	Color
GB-1125601	X-Small	Up to 3.3 cm	White
GB-1125602	Small	3.3 to 4.6 cm	White
GB-1125603	Medium	4.6 to 5.2 cm	White
GB-1125604	Large	5.2 cm to 6.0 cm	White
GB-1125605	X-Large	Over 6.0 cm	White



Measuring point

Use as PIP splint:

L1 = Finger basis fold to distal finger fold

Use as DIP splint:

L2 = PIP joint to distal finger tip

## MobiDig® Extension Finger Splint

Finger splint for extension of PIP joints.  
By carefully bending spring, tension can be increased or decreased.  
Proximal padding positioned in area of palm of hand.  
High wearing comfort by elaborated, wear-resistant padding.  
Linear spring force.

### Indications

Extension of PIP joints


### Sizes

Widths as shown ( 1cm = 0.3937 in).

Reference	Size	L	Finger Width
GB-1125101	Small	6.3 to 7.6	1.5 to 1.7 cm
GB-1125102	Small plus	6.3 to 7.6 cm	1.7 to 1.9 cm
GB-1125103	Medium	7.6 to 8.9 cm	1.7 to 1.9 cm
GB-1125104	Medium plus	7.6 to 8.9 cm	1.9 to 2.1 cm
GB-1125105	Large	8.9 to 10.2 cm	1.9 to 2.1 cm
GB-1125106	Large plus	8.9 to 10.2 cm	2.1 to 2.3 cm



Measuring point:  
L = Distal fold of palm  
of hand to distal finger  
fold



Measuring point:  
Finger width at PIP  
joint




## Capener Flexion Finger Splint

With two lateral spiral springs for anatomically correct movement.  
 Large dorsal pad to spread pressure.  
 Proximal strap for appropriate tension.  
 Progressive spring force.

**Indications:** Flexion of PIP joints

Sizes  
 Widths as shown

Reference	Size	L	Color
GB-1125901	Small	6.3 to 6.9 cm	White
GB-1125902	Medium	6.9 to 7.5 cm	White
GB-1125903	Large	7.5 to 8.2 cm	White
GB-1125904	X-Large	8.2 to 8.5 cm	White



Measuring point:  
 L = Metacarpal knuckle to distal finger fold




## Capener Extension Finger Splint

For extension of PIP joints.  
 With lateral spiral springs for anatomically correct movement.  
 Soft padded.  
 Velcro® closure for fixing.  
 Progressive spring force.

**Indications**  
 Extension of PIP joints

Sizes  
 Widths as shown ( 1cm = 0.3937 in).

Reference	Size	L	Color
GB-1125201	Small	6.3 to 7.0 cm	White
GB-1125202	Medium	7.0 to 7.6 cm	White
GB-1125203	Large	7.6 to 8.0 cm	White




Measuring point  
 L = Distal fold of palm of hand to distal finger fold



## Oval 8® Finger Splint

Oval-8 Finger Splints stabilize and align the small joints of the fingers (DIP and PIP). The versatile Oval-8 splints can be worn six different ways to protect fingers and are ideal to correct deformities caused by arthritis or injury. Conditions include mallet finger, swan-neck deformity, flexible boutonniere deformity, trigger finger, lateral deviation and finger fractures.



### Features:

- Seamless, translucent molded plastic splints are so thin and lightweight you will hardly know they are on your fingers.
- Oval-8 splints do not require strapping or tape to stay in place. Multiple Oval-8 splints can be worn on one hand without interfering with hand function.
- Durable Oval-8 splints can be worn full time or easily slipped on and off as needed. Open design keeps skin dry, yet Oval-8 splints can be worn in water. The broad flat bands distribute pressure evenly for wearing comfort.
- Oval-8 splints can be worn long term by persons with arthritis to align crooked finger joints and may aid in preventing the formation of further deformity.

Oval-8 Sizing Sets are recommended for health care providers to choose and dispense the correct size (s) of Oval-8 Finger Splints for their clients and customers.

Ref.: ALI-51884XX (\*) Size 2 to 15, Full Sizes only


(\*) To order, replace the XX in the reference number by the size, for example: ALI-518843 for size 3 or ALI-5188412 for size 12

## Standard Finger Splints

Used to fix the distal finger joint in the extension position. The hood over the phalanx media should preferably be fixed using the Velcro® strap. The proximal inter-phalangeal joint remains unobstructed. The side walls of the splint are perforated for ventilation.

**Indications:** Mallet Finger, Finger tip and nail bed injuries

**Colors:** Skin tone and clear thermoplastic



Skintone	Clear	Size	Length	Finger Width
PK-102227-00	PK-102228-00	00	4.5 CM	3.5 CM
PK-102227-0	PK-102228-0	0	4.5 CM	4.0 CM
PK-102227-01	PK-102228-01	01	5.5 CM	4.5 CM
PK-102227-02	PK-102228-02	02	5.3 CM	5.0 CM
PK-102227-03	PK-102228-03	03	6.0 CM	5.5 CM
PK-102227-04	PK-102228-04	04	6.2 CM	6.0 CM
PK-102227-05	PK-102228-05	05	6.0 CM	6.5 CM
PK-102227-055	PK-102228-55	5½	6.0 CM	7.0 CM
PK-102227-06	PK-102228-06	06	6.5 CM	7.5 CM
PK-102227-07	PK-102228-07	07	6.2 CM	8.0 CM

## Split Finger Splints

### Adjustable (longitudinal split)

The adjustable, modified finger splints are split right along the dorsal proximal section. Consequently, when the fixing strap is put into position, the splint fits this area of the finger exactly. As the swelling of the finger subsides, the splint may be adjusted easily.

#### Indications:

- Mallet Finger,
- Finger tip
- Nail bed injuries

**Colors:** Skintone

**Sold by the unit**

Reference	Size	Length	Finger Width
PK-102232-00	00	4.5 CM	4.5 CM
PK-102232-0	0	4.5 CM	4.5 CM
PK-102232-01	01	5.5 CM	5.0 CM
PK-102232-02	02	5.3 CM	5.0 CM
PK-102232-03	03	6.0 CM	5.5 CM
PK-102232-04	04	6.2 CM	6.0 CM
PK-102232-05	05	6.0 CM	5.5 CM
PK-102232-055	5½	6.0 CM	7.0 CM
PK-102232-06	06	6.5 CM	7.5 CM
PK-102232-07	07	6.2 CM	8.0 CM

**Ref.: PK-102232-90**

Complete Assortment: 4 pieces of each size above with 20 small Velcro® Straps and 20 large Velcro® straps.



## Boutonnière Finger Splints

This splint has been specially designed for fingers with "buttonhole deformity" (Boutonnière Deformity).

**Indications:** Boutonnière Deformity

**Colors:** Skin tone or Clear thermoplastic

The clear thermoplastic finger splints are made of acrylic resin. The may therefore be reshaped over a gentle flame (spirit or Bunsen burner). This is particularly advantageous for patients whose anatomi-



Color		Size	Length	
Skintone	Clear			
PK-102218-01	PK-102219-01	01	5.0 CM	4.5 CM
PK-102218-02	PK-102219-02	02	5.0 CM	5.0 CM
PK-102218-03	PK-102219-03	03	5.0 CM	5.5 CM
PK-102218-04	PK-102219-04	04	5.0 CM	6.5 CM
PK-102218-05	PK-102219-05	05	5.4 CM	7.0 CM
PK-102218-06	PK-102219-06	06	5.4 CM	7.5 CM

## Digisoft® Finger Orthosis

Finger splint for fixation of 1 or 2 fingers.  
 Pad with soft, ventilating material.  
 Alu-splint can be remodelled.  
 3 Velcro® closures for optimal fixation, center Velcro® closure can be placed in variable position.

### Indications

After injuries and operations, immobilization after plaster cast.  
 Redressement.

Sizes - Widths as shown ( 1cm = 0.3937 in).



Measuring point:  
 L = Beginning of palm to end of finger

Reference	Size	L	Color
GB-1127701	1	Up to 14 cm	Blue
GB-1127702	2	Over 14 cm	Blue

## SellaFix® K Thumbstall

Short version in the carpal joint area for free mobility of wrist. Can be re-worked individually with hair-dryer and scissors. Perforation provides for venting. Adjustable by Velcro® closure. Terry padding strips can be positioned individually. Inner face completely covered with fabric.

### Indications

- For immobilization of
- Metacarpal phalangeal joint with full mobility of wrist.
  - Strains and sprains.

### Sizes

Width of hand at metacarpal knuckle joint as shown ( 1cm = 0.3937 in).



Measuring point  
 H = Circumference of wrist proximal of styloideus

Reference(*)	Size	Circumference	Color
GB-1120801X	X-Small	Up to 15 cm	Beige
GB-1120802X	Small	15-17 cm	Beige
GB-1120803X	Medium	17-19 cm	Beige
GB-1120804X	Large	19-21 cm	Beige
GB-1120805X	X-Large	21-23 cm	Beige

(\*) To order, replace X in the reference number by L for Left, or R for right

## SellaFix® N MetaCarpus Splint

Thumb saddle joint orthosis made of thermoplastic material. Easy correction with hair-dryer and scissors. This allows individual fitting. Fine venting perforation. Washable. Terry padding can be individually positioned. Inner face completely covered with fabric.

### Indications

For immobilization of

- Metacarpal joint and
- Thumb saddle joint.


Ulnar ligament lesion (ski thumb), partial immobilization, rhizo-arthritis.



### Sizing

Width of hand at metacarpal knuckle joint as shown ( 1cm = 0.3937 in).

Reference(*)	Size	H	Color
GB-1120401X	X-Small	Up to 15 cm	Beige
GB-1120402X	Small	15-17 cm	Beige
GB-1120403X	Medium	17-19 cm	Beige
GB-1120404X	Large	19-21 cm	Beige
GB-1120405X	X-Large	21-23 cm	Beige



*Measuring point*

*H = Circumference of wrist proximal of styloideus*

(\*) To order, replace X in the reference number by L for Left, or R for right

## SellaFix® P Thumb

Plus extended thumb tube and Velcro® closure for immobilization of thumb joint, can be readjusted individually with hair-dryer and scissors. Perforation provides for venting. Adjustable by circular Velcro® closure, with terry padding strips for individual positioning. Inner face completely covered with fabric.


### Indications

For immobilization of

- Metacarpal phalangeal joint,
- Thumb saddle joint and
- Interphalangeal joint.
- Tendopathies, fractures.

### Sizing

Circumference of wrist as shown ( 1cm = 0.3937 in).



Measuring point  
H = Circumference of wrist proximal of styloideus

Reference	Size	Circumference	Color
GB-1129001X	X-Small	Up to 15 cm	Beige
GB-1129002X	Small	15-17 cm	Beige
GB-1129003X	Medium	17-19 cm	Beige
GB-1129004X	Large	19-21 cm	Beige
GB-1129005X	X-Large	21-23 cm	Beige

(\*) To order, replace X in the reference number by L for Left, or R for right

## SellaDur Thumb Orthosis

Thumb saddle joint orthosis made of special tough Duroplast. With perforation, therefore very comfortable to wear. In the critical thumb-forefinger-area special high formed edge to avoid pressure spots. Pads to be individually positioned.

### Indications

For immobilization of

- Metacarpal phalangeal joint and
- Thumb saddle joint.


Ulnar ligament lesions (ski thumb), partial immobilization, rhizoarthrosis.

### Sizing

Wrist circumference as shown ( 1cm = 0.3937 in).



Reference(*)	Size	Circumference	Side	Color
GB-1054001X	Small	15-17 cm	Right	Blue
GB-1054002X	Medium	17-19 cm	Right	Blue
GB-1054003X	Large	19-21 cm	Right	Blue



## SellaDur FLOC

Thumb saddle joint orthosis with excellent fit. Perforation provides for venting. In the critical thumb-forefinger-area special high formed edge to avoid pressure spots. Additional flocking with textiles provides optimal wearing comfort. Washable.

### Indications

For immobilization of

- Metacarpal phalangeal joint and
- Thumb saddle joint.


Ulnar ligament lesions (ski thumb), partial immobilization, rhizoarthrosis.

### Sizing

Wrist circumference as shown ( 1cm = 0.3937 in).



Reference(*)	Size	Circumference	Color
GB-1056001X	Small	15-17 cm	Grey
GB-1056002X	Medium	17-19 cm	Grey
GB-1056003X	Large	19-21 cm	Grey



(\*) To order, replace X in the reference number by L for Left, or R for right



## Soft Thumb Splint, Short

Comfortable soft thumb splint. Dry wearing feeling due to naturally breathing and humidity draining material. Integrated anatomically pre-shaped splint for partial immobilization of thumb. Free mobility of the fingers.

Length: Approx. 11 cm


### Indications

After thumb fractures, arthritis, arthrosis of metacarpal phalangeal joint, ulnar ligament lesions (ski thumb).

### Sizing

Widths as shown ( 1cm = 0.3937 in).

Reference	Size	H	Color
GB-1127101	Small	15 to 17 cm	Blue/Black
GB-1127102	Medium	17 to 19 cm	Blue/Black
GB-1127103	Large	19 to 21 cm	Blue/Black



Measuring point  
H = Circumference of wrist proximal of styloideus

## Soft Thumb Splint, Long

Comfortable soft thumb splint. Dry wearing feeling due to naturally breathing and humidity draining material. Integrated anatomically pre-shaped splint for partial immobilization of thumb. Free mobility of fingers.

Extended splint up to wrist area.

Length: approx. 14 cm


### Indications

After thumb fractures, arthritis, arthrosis of thumb saddle joint and/or metacarpal phalangeal joint, ulnar/radial ligament injury of thumb. Prevents false movements, late stage after sinew reconstruction, tendopathy.

### Sizing

Widths as shown ( 1cm = 0.3937 in).

Reference	Size	H	Color
GB-1127201	Small	15 to 17 cm	Blue/Black
GB-1127202	Medium	17 to 19 cm	Blue/Black
GB-1127203	Large	19 to 21 cm	Blue/Black



Measuring point  
H = Circumference of wrist proximal of styloideus