


PLASTER TRAP, VACUUM PRESS, VACUUM FORMING

Q-TRAP (Disposable Plaster Trap)

A unique solution to prevent particulate matter from making its way into the plumbing system. This simple, economical device replaces your existing plaster trap using simple hand tools. Simply remove P-trap from the tailstock in your sink. Attach the hose with the clamps to the sink tailstock, clamping tightly. Take the other hose and insert the tube into the wall drain opening. This should be a snug fit. The large 5-gallon device will stop-up the sink when you need to change the trap. The trap is completely disposable. There is no need to open the trap up: Simply dispose of it and replace with a complete new trap. Tubes measure 1 1/2" in diameter.


Ref.: PEL-QUA591500 5 Gallon
 Ref.: PEL-QUA591501 3.5 Gallon



PLASTER TRAP

- Easy to install and maintain
- Durable non-corrosive all metal construction
- Top mounted drain valve
- Plaster accumulates in a convenient lift-out disposable plastic bag
- Dimensions: 10 1/2" x 9 1/4"
- Capacity: 2 gallons
- Inlet / Outlet openings: 1 1/2"

Ref.: PEL-PLASTERTRAP
 Ref.: PEL-PLASTERTRAPLINERS



VACUUM PRESS 183000

Universal table thermoforming unit for the orthopedic field.

Dimensions:

Outside measurements (LxWxD):
 650MM x 360MM x 150MM

Inside measurements (LxWxD):
 400MM x 180MM x 65MM

Operating Pressure:

Maximum: 8 Bars
 Minimum: 4 Bars

Weight: 9 KG

Ref.: EK-183000



PLASTER TRAP, VACUUM PRESS, VACUUM FORMING

VACUUM PRESS

Quick Vacuum Pull. Table or bench mount.
Complete with hoses, release vacuum valve, and 3/8 hose connectors.

Dimensions: 22" x 23"

Vacuum pump ordered separately.



Ref.: PEL-ROUTMAN-VPT Table Mount
Ref.: PEL-ROUTMAN-VPB Bench Mount

RUBBER SHEETING

For replacement of vacuum forming bladders.

Width: 36"
Sold by the foot



Ref.: PEL-211020 1/8" thick
Ref.: PEL-211008 3/32" thick
Ref.: PEL-211018 1/8" thick

MOBILE VACUUM UNIT

- Sturdy steel frame with casters and wood top.
- Multifunction table-top for pip-mounted platens, height adjustable, vacuum forming directly on table-top and pip-holding fixture for horizontal vacuum forming.
- Master cu-off valve and 3-way valve
 1. Neutral (off) position
 2. Hi vac line-full vacuum direct from vacuum source
 3. Low vac line-controlled by regulator
- On/Off foot operated valve
- Filter and water trap (wet plaster models release 4-6 oz. of water which could contaminate regulator, gauges and pump)
- Six-gallon surge tanks

Vacuum pump ordered separately.

Ref.: PEL-227305 Mobile Vacuum Unit
Ref.: PEL-227005 Vacuum Pump, 110V, 1/4HP



AUTOMATIC VACUUM SYSTEM

- Durable, reliable, dependable.
- 3/4HP oil-less pump.
- 30-gallon tank. Automatic on/off.
- 16-25 inches of mercury

Vacuum pump ordered separately.

Ref.: PEL-ROUTMAN-AVS Single tank
Ref.: PEL-ROUTMAN-AVT Twin tank